


# PRINTER



GLP / GMP Conform documentation of your results.

It provides easily readable records on normal paper which last for years.

It's universal interface allows you to attach "CITIZEN" Balances & Scales via RS232C Interface.

Complete documentation of all data including calibration, weighing results, correct date and time information even after power failure.

Weighing results can be statistically evaluated. Mean value and standard deviation as well as the total weight are automatically printed out due to built-in statistics and totalization function.

## Technical Specification :

Model	CPR02 Date, time, statistic and totalization functions
Print	Dot Matrix Printer
No of characters per line	24 Character
Printing Speed	Faster than 1 line per second
Paper	Normal Paper Roll, Wood Free Paper, Width: 57.5 ± 0.5mm, Core Diameter : 8mm Min.
Power Supply	AC Input 100V - 300V & 50Hz DC Output 13V & 1.5 Amps
Interface	RS 232 C
Transmission parameters	Speed of transmission : 1200, 2400, 4800, 9600, 19200 baud. Data Bit - 8, Stop Bit 1, Parity - None, Space, Odd, Even
Working Conditions	Working temperature : 5°C to 50°C, relative humidity: 10% - 80% (no condensation)
Masses and dimensions (mm)	88.5 (H) x 135 (W) x 210 (L)
Components	9 pin Rs232 Cable, 13V DC Adapter, Ink Ribbon Paper Roll, user manual and warranty

### CITIZEN SCALE (I) PVT. LTD.

Unit No. E2, Plot No. 11, WICEL Estate, Opp. Seepz Gate No.1, Andheri (E), Mumbai 400 093. (INDIA)  
Tel. : +91-22-4243 7700 • Fax : +91-22-4243 7800 • E:mail : sales@citizenscales.com • Web : www.citizenscales.com