

DIGITAL ULTRASONIC CLEANER / BATH



Citizen®
ISO 9001 : 2000 CERTIFIED

10 Ltr. With Heater

YJ 5200 DT (S.S.)



with
Steam Fresh

Features :

- Bigger tank with Special basket for cleaning.
- M.S. & S.S. Housing available for rugged conditions.
- Digital LED / LCD digital for smooth working.
- Ideal for any kind of difficult conditions.
- Heater & Timer Function available.
- Water spout for drainage.
- Strong transducer for years of operations.
- S.S. 304 inside Cleaning Tank for all Models.



Waterspout



Special basket for cleaning



Adjust Time and Temp.



Family



Restaurant



Repair Shop



Pharma



Jewellery Factory



Chemical Lab



Dental Clinic

Technical data:

MODEL	YJ 5200 - DT
Volume (Ltr.)	10Ltr.
Voltage	220V / 50Hz - 110V / 60Hz
Ultrasonic Power (W)	240W
Heat Power	225W
Adj. Power (%)	20 - 80
Temp. Range (°C)	1 - 60
Frequency (KHz)	40 Khz
Tank Size (mm)	300 x 240 x 150
Unit Size (mm)	420 x 350 x 420
G.W.	10kg
Indication	LED Digital Indication

WORLDWIDE MARKETING:

CITIZEN SCALE (I) PVT. LTD.

Unit No. E2, Plot No. 11, WICEL Estate, Opp. Seepz Gate No.1, Andheri (E), Mumbai 400 093. (INDIA)
Tel. : +91-22-4243 7700 • Fax : +91-22-4243 7800 • E:mail : sales@citizenscales.com • Web : www.citizenscales.com